

This package insert must be read carefully prior to use.

Total protein assay kit for blood tests  
(Classification No.:30181001)

Biuret reagent  
Protein in the sample  $\longrightarrow$  Complex  
(Red-purple color)

## Clinimate TP

### General Precautions \*\*

1. This product is for in vitro diagnostic use, and must not be used for any other purposes.
2. Clinicians should make a comprehensive clinical decision based on assay results in conjunction with clinical symptoms and other examination results.
3. For the effects of an administered drug on the measured value, carefully read the Precautions for Use in the package insert of the drug, especially the section about the effects on laboratory test results. Please also read carefully the “2. Interfering Substances,” in the “Procedural Precautions” section, as well as “2. Precautions for Assessment” in the “Assessment of Assay Results” section, of this package insert.
4. This product should be used only as directed in this package insert. Reliability of results cannot be guaranteed if there are any deviations from the instructions in this package insert.
5. If the reagent accidentally comes in contact with eyes and/or mouth, rinse immediately with ample water as first aid, and consult the doctor if required.
6. Carefully read the operating instructions for each type of automated analyzers prior to using this product. Parameters for each type of analyzers are available, and can be requested from SEKISUI MEDICAL CO., LTD. if required.
7. Perform a quality control test prior to assay to ensure accuracy.

### Description (Kit Components) \*\*

Component: Ingredients

TP Coloring Solution:

Copper sulfate (II) pentahydrate

### Intended Use

#### Measurement of total protein in serum or plasma

Total protein is the general term for all proteins with various functions. Changes in total protein levels reflect change of physiological function and damage or disease affecting organs/tissues of the body.

### Assay Principle

#### 1. Assay Principle

Protein in the samples reacts with the Biuret reagent to form a red-purple complex. The amount of total protein is determined by measuring the absorbance of this red-purple complex.

### 2. Features

- 1) Single liquid method without any need for reagent preparation.
- 2) The measurement is hardly affected by turbidity of the solution.

### Procedural Precautions \*\*

#### 1. Properties of Samples and Sampling Methods

- 1) Samples  
Serum and plasma (heparin plasma, EDTA plasma, citrate plasma and NaF-EDTA plasma) may be used.
- 2) Storage of samples<sup>3)</sup>  
If the isolated serum or plasma sample cannot be tested on the same day, specimens should be stored as follows:  
2–10°C: for tests within 1 month  
≤ -20°C: for tests after more than 1 month  
Bring samples to room temperature (15–30°C) before use.

#### 2. Interfering substances

Icteric and hemolyzed blood samples cause negative and positive measurement errors, respectively. Perform re-measurement by another test method if the results obtained are extremely low or high compared with those of other related variables.

#### 3. Others

- 1) Always use Serum Multicalibrator (SEKISUI), Seronorm Multicalibrator, Seronorm Human or Anaserum ALB/TP Standard Solution for calibration.
- 2) Analysis range  
If the concentration of sample exceeds assay range, dilute the sample with saline and repeat the measurement.

### Dosage/Administration (Assay Procedure) \*\*

#### 1. Preparation of reagents

Reagent: TP Coloring Solution is ready to use.

#### 2. Assay Procedure

This product is compatible with various types of automated analyzers. An example of the assay procedure is indicated below.

Sample 5 μL	+	Reagent 350 μL	$\xrightarrow[7.5 \text{ min}]{37^\circ\text{C}}$	Measurement (Absorbance <sup>**</sup> )
				$\longrightarrow$
				Calculation of concentration

<sup>\*\*</sup>Absorbance: The difference in absorbance between 700 nm and 546 nm  
Calibration material: Serum Multicalibrator (SEKISUI), Seronorm Multicalibrator, Seronorm Human or Anaserum ALB/TP Standard Solution (Manufacturer's assigned value)

Reagent blank: Purified water or saline

### Assessment of Assay Results \*\*

#### 1. Reference standard range <sup>2)</sup>

6.6–8.1 g/dL (Within the JCCLS common standard)

#### 2. Precautions for Assessment

There may be reactions or interfering reactions with non-target substances. If assay results appear to be unreliable, repeat the measurement (if necessary, after dilution) or try another analytical methods.

### Performance \*

#### 1. Sensitivity

- 1) Reagent blank: absorbance being equal to or lower than 0.05
- 2) Sensitivity: The absorbance is 0.21–0.25 per 6 g/dL of albumin.

#### 2. Accuracy: 90–110% of the expected assay value

#### 3. Within-run reproducibility:

Coefficient of variation  $\leq 3\%$   
(Test methods used for 1. –3. are in-house methods.)

#### 4. Measurement Range<sup>5)</sup>: (On Hitachi 7170S automated analyzer)

0.1–13.0 g/dL

#### 5. Correlation<sup>5)</sup>

- 1) Serum N=60  $r=0.972$   $y=1.08x - 0.61$   
Control method: Approved in vitro diagnostic (Biuret test)
- 2) Plasma N=150  $r=0.999$   $y=1.02x - 0.11$   
Control method: Approved in vitro diagnostic (Biuret test)

#### 6. Standard Material

SRM927 (NIST)

### Precautions for Use or Handling \*

#### 1. Precautions for Handling (to Ensure Safety)

- 1) All samples used in the test should be handled as a material possibly infected with HIV, HBV, HCV, or other viruses. To prevent infection, use disposable gloves and avoid mouth pipetting during the test.
- 2) The TP Coloring Solution contains sodium hydroxide. If it is accidentally ingested or comes into contact with the eyes or skin, immediately rinse with water and seek medical treatment, if necessary.

#### 2. Precautions for use

- 1) This product should be stored as directed, avoid freezing. Freezing can cause deterioration of the reagents, leading to inaccurate results. Therefore, do not use the product if it has been previously frozen.
- 2) Do not use expired reagents. Use of such reagents cannot guarantee the reliability of measurement values.
- 3) Do not replenish the reagents.
- 4) Do not perform the assay under direct sunlight.

#### 3. Precautions for Disposal

- 1) Before disposal, used samples and their containers must be immersed in sodium hypochlorite solution at a concentration of greater than 0.1% for longer than 1 hour or autoclaved at 121°C for 20 minutes.

- 2) To prevent infections from spilled samples or solutions containing samples, wipe the spilled area thoroughly with disinfectants such as sodium hypochlorite solution at a concentration of greater than 0.1%.
- 3) The reagents and treated samples should be discarded as medical waste or industrial waste according to the waste disposal regulations.
- 4) The reagents should be disposed of in accordance with the Water Pollution Control act or related regulations.

#### 4. Other precautions

Do not use the containers for other purposes.

### Storage and Shelf Life \*

1. Storage temperature: room temperature
2. Shelf life: 2 years from the date of manufacture (The expiration date is printed on the outer package.)

### Packaging

Name		Package
Clinimate TP	TP Coloring Solution	4 × 100 mL

Constituent reagents are available in other configurations. For further details please contact SEKISUI MEDICAL CO., LTD.

### References \*\*

- 1) Kitamura M. et al. Practical Clinical Chemistry. Tokyo: Ishiyaku Shuppan; 1974. Japanese.
- 2) Kanai M, editor. Kanai's manual of clinical laboratory medicine. 35th ed. Tokyo: Kanehara Shuppan; 2020. Japanese.
- 3) Sasaki M. et al. Sampling of constituents of the human body. Tokyo: Kodansha; 1972. Japanese.
- 4) Nakahara M. et al.: Jpn J Med Tech, 45, 1106, 1996. Japanese.
- 5) SEKISUI MEDICAL CO., LTD. In house data. Japanese.

### Contact

SEKISUI MEDICAL CO., LTD.  
international@sekisui.com

### Manufacturer \*

**SEKISUI MEDICAL CO., LTD.**

**1-3, Nihonbashi 2-chome, Chuo-ku, Tokyo, Japan**

"Clinimate" and "Anaserum" are trademarks or registered trademarks of SEKISUI MEDICAL CO., LTD. in Japan and/or other countries.

All other brands, product names and service names are trademarks or registered trademarks of their respective companies.