

For in vitro diagnostic use only.

Revised: January 2017

HbA1c Pretreatment Solution for NORUDIA N HbA1c

Instructions for Use

This product is used to pretreat blood cell samples for the quantitative determination of Hemoglobin A1c (HbA1c) in whole blood using NORUDIA N HbA1c.

Preparation

1. This product is stable for 4 weeks at room temperature after opening.
2. Extract 25 μ L of blood cell layer from whole blood (EDTA, NaF-EDTA, Heparin) that has been stored overnight or centrifuged (800 g for 5 minutes) and add it to 500 μ L of this Pretreatment Solution. Mix to prepare the pretreated sample.
3. The pretreated assay sample should be tested within 8 hours at room temperature or within 24 hours when refrigerated at 2 - 10°C. Such stored samples should be brought to room temperature (15–30°C) before use.
4. Mix the sample well before the assay if stored for more than 2 hours.

Precautions

1. Use this Pretreatment Solution for NORUDIA N HbA1c. Do not use it for any other purpose.
2. This product contains Proclin 300 which may irritate skin. If the reagent accidentally comes into contact with the skin, clothing or enters the mouth or the eyes, flush with water immediately, and seek medical attention accordingly.
3. Disposal of the remaining reagents, samples and accessories should be carried out according to the regulations for waste products and handled separately from medical and industrial waste.

Packaging

6 \times 50 mL

Storage

2–30°C

Expiration Date

Printed on the outer package

SEKISUI MEDICAL CO., LTD.

1-3, Nihonbashi 2-chome, Chuo-ku, Tokyo, Japan