

# パッケージ試験 測定項目

## 探索スクリーニングパッケージと副次的薬理試験パッケージ

P1,2,3:探索スクリーニング、副次:副次的薬理試験

	測定系	サブタイプ	レセプター	P1	P2	P3	副次
1	Adenosine	A1	Human				
2		A2a	Human				
3		A3	Human				
4		Transporter	Human				
5	Adrenergic	α1A	Rat				
6		α1B	Rat				
7		α2A	Human				
8		α2B	Human				
9		α2C	Human				
10		β1	Human				
11	β2	Human					
12	Angiotensin	AT1	Human				
13		AT2	Human				
14	Bradykinin	B1	Human				
15		B2	Human				
16	Bombesin	Non-Selective	Rat				
17	Ca channel	Type L (Benzothiazepine)	Rat				
18		Type L (Dihydropyridine)	Rat				
19		Type L (Phenylalkylamine)	Rat				
20		Type N	Rat				
21	Cannabinoid	CB1	Human				
22		CB2	Human				
23	CCK	A	Human				
24		B	Human				
25	CRF	CRF1	Human				
26	Dopamine	D1	Human				
27		D2 short	Human				
28		D3	Human				
29		D4.2	Human				
30		D5	Human				
31		Transporter	Human				
32	Estrogen		Rat				
33	Endothelin	ETA	Human				
34		ETB	Human				
35	GABA	GABA A (Agonist Site)	Rat				
36		GABA A (BZ Central)	Rat				
37		Chloride Channel	Rat				
38		GABA B	Rat				
39		Transporter	Rat				
40	Glucocorticoid		Human				
41	Glutamate	AMPA	Rat				
42		Kainate	Rat				
43		NMDA (Agonist site)	Rat				
44		NMDA (Glycine site)	Rat				
45		NMDA (Phencyclidine site)	Rat				
46		NMDA (Polyamine site)	Rat				

	測定系	サブタイプ	レセプター	P1	P2	P3	副次
47	Glycine	Strychnine sensitive	Rat				
48	Histamine	H1	Human				
49		H2	Human				
50		H3	Human				
51	Imidazoline	Central	Rat				
52	IP3		Rat				
53	K channel	KA	Rat				
54		KATP	Rat				
55		Skca	Rat				
56	Leukotriene	B4	Guinea pig				
57		D4	Guinea pig				
58	Melatonin	MT1	Human				
59	Monoamine	Transporter	Rabbit				
60	Muscarinic	M1	Human				
61		M2	Human				
62		M3	Human				
63		M4	Human				
64		M5	Human				
65	Na channel		Rat				
66	Neurokinin	NK1	Human				
67		NK2	Human				
68		NK3	Human				
69	Neuropeptide	NPY1	Human				
70		NPY2	Human				
71	Neurotensin	NT1	Human				
72	Norepinephrine	Transporter	Human				
73	Nicotinic	Ni	Human				
74	Opiate	δ	Human				
75		κ	Human				
76		μ	Human				
77		ORL1	Human				
78	PAF		Rabbit				
79	Prostanoid	EP2	Human				
80	Serotonin	5HT1A	Human				
81		5HT2A	Human				
82		5HT3	Human				
83		Transporter	Human				
84	Sigma	σ1	Guinea pig				
85		σ2	Guinea pig				
86	Testosterone		Human				
87	Vasopressin	V1	Rat				
88		V1B	Human				
89		V2	Human				
90	VIP	VIP	Human				

測定項目数 20 60 90 60