

◆Chemicals

品目コード	メーカーコード	品名		包装	貯法
497445	CYP500	4-hydroxydiclofenac, 0.1 mg	CYP2C9代謝物	0.1 mg	-20°C
497452	CYP500-5	4'-hydroxydiclofenac, 5 mg	CYP2C9代謝物	5 mg	-20°C
497469	CYP501	6 α -hydroxypaclitaxel, 0.1 mg	CYP2C8代謝物	0.1 mg	-20°C
497476	CYP501-5	6 α -hydroxypaclitaxel, 5 mg	CYP2C8代謝物	5 mg	-20°C
497513	CYP505	Tolbutamide, 100 mg	CYP2C基質	100 mg	室温
497520	CYP506	Hydroxytolbutamide, 10 mg	CYP2C代謝物	10 mg	室温
497537	CYP507	Carboxytolbutamide, 10 mg	CYP2C代謝物	10 mg	室温
497544	CYP508	Carboxytolbutamide methyl ester, 10 mg		10 mg	室温
497551	CYP509	Tolbutamide kit (CYP505 - CYP508)		1 kit	室温
497568	CYP510	7-ethoxyresorufin, 10 mg	CYP1A/1B基質	10 mg	-20°C
497575	CYP511	Resorufin, 100 mg	CYP1A/1B 代謝物	100 mg	-20°C
497582	CYP512	7-benzyloxyquinoline, 10 mg	CYP3A4基質	10 mg	-20°C
497599	CYP513	7-hydroxyquinoline, 10 mg	CYP3A4代謝物	10 mg	室温
497612	CYP515-10	MAMC, 10 mg	CYP2D6基質	10 mg	-20°C
497629	CYP515-25	MAMC, 25 mg	CYP2D6基質	25 mg	-20°C
497636	CYP516	HAMC, 10 mg	CYP2D6代謝物	10 mg	-20°C
497643	CYP517-10	7-methoxy-4-(trifluoromethyl)-coumarin (MFC), 10 mg	CYP2C9基質	10 mg	-20°C
497650	CYP517-25	7-methoxy-4-(trifluoromethyl)-coumarin (MFC), 25 mg	CYP2C9基質	25 mg	-20°C
497667	CYP518-10	7-benzyloxy-4-(trifluoromethyl)-coumarin (BFC), 10 mg	CYP3A4基質	10 mg	-20°C
497674	CYP518-25	7-benzyloxy-4-(trifluoromethyl)-coumarin (BFC), 25 mg	CYP3A4基質	25 mg	-20°C
497681	CYP519	7-hydroxy-4-(trifluoromethyl)-coumarin (HFC), 10 mg	CYP3A4/2C9 代謝物	10 mg	-20°C
497698	CYP520	Omeprazole, 1 g	CYP2C19基質	1 g	2~8°C
497704	CYP521-10	3-O-methylfluorescein, 10 mg	CYP2C19基質	10 mg	-20°C
497711	CYP521-25	3-O-methylfluorescein, 25 mg	CYP2C19基質	25 mg	-20°C
497728	CYP522-10	Tienilic Acid, 10 mg	CYP2C9阻害剤	10 mg	-20°C
497735	CYP522-25	Tienilic Acid, 25 mg	CYP2C9阻害剤	25 mg	-20°C
497742	CYP523-25	Paclitaxel, 25 mg	CYP2C8基質	25 mg	-20°C
497759	CYP523-100	Paclitaxel, 100 mg	CYP2C8基質	100 mg	-20°C
497766	CYP524-100	Ketoconazole, 100 mg	CYP3A4阻害剤	100 mg	2~8°C
497773	CYP524-1000	Ketoconazole, 1 g	CYP3A4阻害剤	1 g	2~8°C
497780	CYP525	Diclofenac, 10 g	CYP2C9基質	10 g	室温
497797	CYP526	Amodiaquine dihydrochloride dihydrate, 5 g	CYP2C8基質	5 g	室温
497803	CYP527-10	Desethylamodiaquine, 10 mg	CYP2C8代謝物	10 mg	室温
497810	CYP527-25	Desethylamodiaquine, 25 mg	CYP2C8代謝物	25 mg	室温
497827	CYP528	Terfenadine, 5 g	CYP2J2/3A4 阻害剤	5 g	2~8°C
497834	CYP529-10	R-(-)-N-benzylphenobarbital, 10 mg	CYP2C19阻害剤	10 mg	-20°C
497841	CYP529-25	R-(-)-N-benzylphenobarbital, 25 mg	CYP2C19阻害剤	25 mg	-20°C
497858	CYP530-10	Montelukast sodium, 10 mg	CYP2C8阻害剤	10 mg	室温
497865	CYP530-25	Montelukast sodium, 25 mg	CYP2C8阻害剤	25 mg	室温
497872	CYP531-10	Diethoxyfluorescein (DEF), 10 mg	CYP2C8/3A4 基質	10 mg	室温
497889	CYP531-25	Diethoxyfluorescein (DEF), 25 mg	CYP2C8/3A4 基質	25 mg	室温
497896	CYP532	S-mephenytoin, 5 mg	CYP2B6/2C19 基質	5 mg	2~8°C
497902	CYP532-25	S-mephenytoin, 25 mg	CYP2B6/2C19 基質	25 mg	2~8°C
497919	CYP533	(\pm)-4'-hydroxymephenytoin, 5 mg	CYP2B6/2C19 代謝物	5 mg	2~8°C
522536	CYP534	1'-hydroxymidazolam, 10 mg	CYP3A代謝物	10 mg	2~8°C
522543	CYP535	4-hydroxymidazolam, 10 mg	CYP3A代謝物	10 mg	2~8°C